

# The Multi-Task Genie™ Family



**HIGH SPEED MIXER**



**THE ORIGINAL**



**INTEGRATED TIMER**



**DIGITAL CONTROL**

Vortex-Genie® 1	Vortex-Genie® 2	Vortex-Genie® 2T	Digital Vortex-Genie® 2
<ul style="list-style-type: none"> <li>Violent single high-speed vortexing</li> <li>Easy-to-use Touch On operation</li> <li>Large, 6mm orbit for aggressive vortexing of difficult samples</li> <li>May be used with Pop-off Cup (shown) or optional 3-inch Platform</li> <li>May be used in cold rooms and incubators</li> </ul>	<ul style="list-style-type: none"> <li>OUR MOST POPULAR MODEL</li> <li>Variable speed for precise mixing from gentle to vigorous</li> <li>Broad range of attachments for most tubes, plates and other containers</li> <li>Hands-Free or Touch On control</li> <li>May be used in cold rooms and incubators</li> </ul>	<ul style="list-style-type: none"> <li>Same great features as the Vortex-Genie 2 with an integrated timer</li> <li>Timer function:                             <ul style="list-style-type: none"> <li>- Touch On: 1-60 sec.</li> <li>- Hands-Free: 1-60 min. or continuous</li> </ul> </li> <li>Unit automatically stops at the end of user-programmed time</li> <li>May be used in cold rooms and incubators</li> </ul>	<ul style="list-style-type: none"> <li>Same great features as the Vortex-Genie 2 with digital control and display of speed and time</li> <li>Digital display provides accuracy, reproducibility and repeatability</li> <li>Timer function:                             <ul style="list-style-type: none"> <li>- Touch On: 1-99 sec.</li> <li>- Hands-Free: 1-99 min. or continuous</li> </ul> </li> <li>May be used in cold rooms and incubators</li> <li><b>Specially configured unit for Applied Biosystems protocol available</b></li> </ul>

## SPECIFICATIONS

Vortex-Genie Model	Vortex-Genie 1	Vortex-Genie 2	Vortex-Genie 2T	Digital Vortex-Genie 2
Supplied	Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup	3-inch Platform, Pop-off Cup
Speed (RPM)	3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	600-3200 (2700 for 50Hz models)	500-3000 (2850 for 50Hz models)
Timer - Touch On: Hands-Free:	N/A	N/A	1-60 seconds 1-60 minutes or continuous	1-99 seconds 1-99 minutes or continuous
Weight	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)	4 Kg (8.8 lb)
Dimensions (base) (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)
Markings / Warranty	ETL, CE / 2 years	ETL, CE / 2 years	ETL, CE / 2 years	CE / 2 years

## CATALOG NUMBERS

				AB MODELS	
120V (North America)	SI-0136	SI-0236 (Model G560)	SI-T236	SI-A236	SI-A536
230V - No Plug	SI-0146	SI-0246 (Model G560E)	SI-T246	SI-A246	SI-A546
230V - European Plug	SI-0156	SI-0256 (Model G560E)	SI-T256	SI-A256	SI-A556
230V - British Plug	SI-0166	SI-0266 (Model G560E)	SI-T266	SI-A266	SI-A566
230V - Swiss Plug	SI-0176	SI-0276 (Model G560E)	SI-T276	SI-A276	SI-A576
240V - Australian Plug	SI-0196	SI-0297 (Model G560E)	SI-T296	SI-A296	SI-A596
100V	SI-0186	SI-0286	SI-T286	SI-A286	SI-A586

See back page for additional Vortex-Genie 2, Vortex-Genie 2T and NEW Digital Vortex-Genie 2 accessories.

U.S. Pats. 4,781,487 / 5,707,861



Home of the Multi-Task Genie™ Family

80 Orville Drive, Suite 102  
Bohemia, NY 11716 U.S.A.  
Phone: (631) 567-4700 • Fax: (631) 567-5896  
Toll Free: 888-850-6208

Customer Service: custsvc@scientificindustries.com  
scientificindustries.com

# Multi-Task Genie™ Family Accessories

Accessories for all Vortex-Genie® 2 mixers and Vortex-Genie® Pulse




Fig. 1  
Model H301  
Multiple  
Sample  
Starter Set




Fig. 4  
Model H302  
Large Sample  
Set



Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	13	SI-0564	TurboMix™ Attachment
--	146-6005-00	6-inch Platform	14	SI-0565	Universal Microtube Holder*
2	504-0234-00	Microtube Foam Insert (2)**	15	SI-0511	Large Ampule/Tube Attachment
3	504-0235-00	Microplate Foam Insert (2)**	16	SI-0570	Small Ampule/Tube Attachment
3A	504-0232-00	Deepwell Microplate Foam Insert (2)**	17	SI-H506	Horizontal 50ml Tube Holder (6 tubes)
4	503-0282-00	Model H302 Large Sample Set	18	SI-H512	Horizontal 15ml Tube Holder (12 tubes)
5	504-0039-00	Recessed Platform	19	SI-H524	Horizontal Microtube Holder (24 tubes)
6	503-0278-00	9-16mm Tube Foam Inserts (2)***	20	SI-V203	Vertical High Speed 50ml Tube Holder (3 tubes)
7	503-0279-00	14-29mm Tube Foam Inserts (2)***	21	SI-V206	Vertical High Speed 15ml Tube Holder (6 tubes)
8	503-0280-00	29-37mm Tube Foam Inserts (2)***	22	SI-V209	Vertical High Speed Microtube Holder (9 tubes)
9	503-0277-00	Blank Foam Inserts (2)***	23	SI-V506	Vertical 50ml Tube Holder (6 tubes)
--	503-0283-00	Assorted Foam Inserts (4)***	24	SI-V512	Vertical 15ml Tube Holder (12 tubes)
10	568-0001-00	Elastic Bands (2)***	25	SI-V524	Vertical Microtube Holder (24 tubes)
11	146-3011-00	Pop-off Cup	26	SI-V525	Vertical (6) 50ml, (6) 15ml, (12) Microtube Holder
12	0K-0500-902	3-inch Platform with Rubber Cover	--	0K-0236-408	Bottom Closure with Screws Kit †
--	580-2013-00	Rubber Cover for 3-inch Platform			

\* May be used with TurboMix attachment or Disruptor Genie®

\*\* 6-inch Platform Required

\*\*\* Recessed Platform Required

† May not fit older models. Please inquire.

**NOTE:** Vortex-Genie 1 may only be used with Pop-off Cup and 3-inch Platform

U.S. Pat. 5,707,861



**Home of the Multi-Task Genie™ Family**  
 80 Orville Drive, Suite 102  
 Bohemia, NY 11716 U.S.A.  
 Phone: (631) 567-4700 • Fax: (631) 567-5896  
 Toll Free: 888-850-6208  
 Customer Service: custsvc@scientificindustries.com  
[scientificindustries.com](http://scientificindustries.com)

© 2014 Scientific Industries LTP0049 Rev. H